MVE TEC 3000

Quick Reference User Manual Multilingual

MVE Biological Solutions
3655 Torrington Drive
Ball Ground, GA 30107
+1-678-683-2799
USA www.mvebio.com

Medical Product Services
Bongasse 20
35619 Bräunle, Germany
+49-6942-962073

CE 0459

Intended Use & Indication for Use for Cryogenic Storage and/or Transport


LIFE SCIENCE INTENDED USE STATEMENT

STORAGE ONLY

MVE FREEZERS ARE INTENDED FOR THE INDICATION OF PRESERVING HUMAN BIOLOGICAL PRODUCTS, SAMPLES, AND/OR SPECIMENS (E.G., BLOOD, BLOOD PRODUCTS, CELLS, TISSUES, ETC.) AT CRYOGENIC AND ULTRACOLD TEMPERATURES DURING STORAGE.

STORAGE AND TRANSPORT

MVE DEWARs AND VAPOR SHIPPERS ARE INTENDED FOR THE INDICATION OF PRESERVING HUMAN BIOLOGICAL PRODUCTS, SAMPLES, OR SPECIMENS (E.G., BLOOD, BLOOD PRODUCTS, CELLS, TISSUES, ETC.) AT CRYOGENIC AND ULTRACOLD TEMPERATURES DURING STORAGE AND/OR TRANSPORT.

TEC 3000 Quick Reference User Manual

Table of Contents

Product Identification	
Display / Control Panel.....	2
Bottom Panel / Electrical / Physical Connections.....	3
Dewar Plumbing Connections	4
Adjusting Temperature Alarm Settings.....	5
Adjusting Level and Level Alarm Settings.....	6-7
High Level Alarm Setting	
High Level Setting	
Low Level Setting	
Low Level Alarm Setting	
Adjusting Display and Output Settings.....	8
Calibration of Temperature Probes.....	8
Password and Security Setup.....	9-10
Alarms and Definitions.....	11
Contact Information	11
German table of contents	12
Italian table of contents.....	23
Spanish table of contents	34
French table of contents	45

NOTE: All models have been tested for electrical safety per IEC 60601-1 Class 1 requirements and found to be compliant. Select models are approved and labeled as Class IIa medical devices per the European Union, Medical Device Directive. The US FDA does not classify these models as medical devices.

NOTE: MVE liquid nitrogen freezers should be installed by an authorized MVE Distributor per the TEC 3000 Technical Manual, PN 13289499.

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3. Safety and First Aid

Symbols Used

The following symbols are used in this manual, on the device, and on device packaging:

Symbol	Title	Description
	Operating Instructions	The operating instructions should be considered for additional information when operating this device.
	Serial Number	Unique identifier for the device.
	Model Number	MVE model number for the device.
	Caution	Signifies a CAUTION of a potentially hazardous situation when operating the device that may result in minor to moderate injury or property damage.
	Warning	Signifies a WARNING of a potentially hazardous situation when operating the device that may result in serious injury or property damage.
	Warning; Low Temperature	Indicates low temperature or freezing conditions. Take care to avoid exposure to skin, eyes, and clothing.
	Warning; Asphyxiating Atmosphere	Indicates the potential for an oxygen-depleted atmosphere due to nitrogen vapor. Take care to operate device in a well-ventilated area.
	Warning; Electricity	Indicates a potential electrical hazard. Take care to avoid contact with electricity.
	Warning; Explosive	Indicates a potential explosive hazard. The expansion ratio of liquid nitrogen to gas is 1:700 and can cause explosive conditions if placed into a
	Wear Protective Gloves	Thermal gloves must be worn during indicated procedures.
	Wear a Face Shield	A face shield must be worn during indicated procedures.
	No Pushing	For the MVE 200 series of freezers. Indicates the area of the freezer that should not be pushed due to potential overbalance.
	Temperature Limit	Indicates minimum and maximum temperature limits at which the freezer should be stored or transported.
	Humidity Limit	Indicates minimum and maximum humidity limits at which the freezer should be stored or transported.
	UL Listed Mark	MVE Cryogenic Freezers conform to relevant UL safety standards.
	CE Mark	MVE Cryogenic Freezers are assessed to meet safety, health, and environmental protection requirements for Europe.
	Manufacturer	Indicates manufacturer name and address.

- **WARNING: Do not modify this equipment without authorization of the manufacturer.**

TEC 3000 Quick Reference User Manual

Liquid Nitrogen Safety

Liquid nitrogen (LN2) is used in MVE Cryogenic Freezers as a refrigerant. Understanding potential hazards and following safety precautions is important when handling LN2 and these freezers. Nitrogen is a colorless, odorless, and tasteless gas that makes up approximately 78.1% of the Earth's atmosphere in its gaseous state. LN2 becomes vapor at temperatures greater than -320.8°F (-196°C). In liquid state, nitrogen has a temperature range from -320.4°F to -346°F (-195.8°C to -210°C).

- Nitrogen vapor is a potential asphyxiant as it displaces Oxygen (O₂) in confined spaces. Rapid suffocation can occur without warning in an Oxygen-deficient atmosphere (less than 19.5% O₂). MVE Cryogenic Freezers must be installed and operated in well-ventilated areas.
- DO NOT vent container in confined spaces.
- DO NOT enter confined spaces where excess nitrogen gas may be present.
- If exposure has occurred move to ventilated area or fresh air. If breathing is difficult, supplement oxygen may be required. If not breathing, give artificial respiration. SEEK IMMEDIATE MEDICAL ATTENTION.

- Contact with liquid nitrogen or uninsulated equipment containing nitrogen can result in cold contact burns or tissue damage. Nitrogen vapor can cause damage to skin or eyes.
- In case of frostbite, warm area with warm water not exceeding 105°F (40°C) and SEEK IMMEDIATE MEDICAL ATTENTION.

- Never place LN2 in a sealed container without a pressure relief device. The expansion ratio of liquid nitrogen to gaseous nitrogen is 1 to 700 (1 cubic foot of liquid nitrogen becomes 700 cubic feet of gaseous nitrogen when evaporated).

Recommended protective clothing

- Cryogenic gloves (loose fitting)
- Full-face shield or chemical splash goggles
- Cryogenic apron
- Long sleeve shirt and cuffless pants
- Closed toe shoes (no sandals)

TEC 3000 Quick Reference User Manual

Equipment Usage

Cryogenic containers must be operated in accordance with the manufacturer/supplier instructions. Safety instructions will also be posted on the side of each Dewar. Cryogenic Dewars must be kept in a well-ventilated area protected from weather and away from heat sources. In applications that use a modular liquid cylinder as a source of LN₂, the supply will need to be replenished at regular intervals to ensure proper operation of the freezer. When exchanging liquid cylinders, follow the below procedure:

1. Allow all plumbing components to warm to room temperature before attempting to change supplies.
2. Close all valves associated with the liquid supply cylinder.
3. Relieve pressure in the plumbing assembly by initiating a brief fill by either pressing "Start Fill".
4. Loosen the plumbing connection for the transfer hose at the liquid cylinder.
5. Remove empty liquid cylinder and replace with full liquid cylinder pressurized to 22 - 35 psig (1.52 - 2.41 bar).
6. Attach the transfer hose to the plumbing connection on the liquid cylinder. Ensure that the hose is connected to the connection labeled "LIQUID".
7. Tighten the transfer hose plumbing connection at the liquid cylinder.
8. Open the liquid supply valve on the liquid cylinder.
9. Inspect plumbing for audible and visual leaks. Repair if necessary.
10. Manually initiate a fill to verify proper operation.

Recommended First Aid

Every site that stores and uses LN₂ should have an appropriate Material Safety Data Sheet (MSDS) present. The MSDS may be obtained from the manufacturer/distributor. The MSDS will specify the symptoms of overexposure and first aid to be used. Here is a typical summary.

- If symptoms of asphyxia such as headache, drowsiness, dizziness, excitation, excess salivation, vomiting, or unconsciousness are observed, remove to fresh air. If breathing has stopped, give artificial respiration. **CALL A PHYSICIAN IMMEDIATELY.** If breathing is difficult, supplemental oxygen may be required.
- If exposure to cryogenic liquids or cold vapor occurs, restore tissue to normal, body temperature (37°C) as rapidly as possible, and then protect the injured tissue from further damage and infection. Rapid warming of the affected areas is best achieved by bathing it in warm water. The temperature of the water used should not exceed 40°C. Under no circumstances should the frozen part be rubbed either before or after warming. If the eyes are involved, flush them thoroughly with warm water for at least 15 minutes. In case of massive exposure, remove clothing while showering with warm water. The patient should not drink alcohol or smoke. **CALL A PHYSICIAN IMMEDIATELY.**

TEC 3000 Quick Reference User Manual


Table 1: Front Panel Identification

NOTE: The recommended LN2 supply pressure is 22-35 psig (1.52 - 2.41 bar).

1	Front Panel	The front panel is the user interface for the TEC 3000. All displays and controls are located on the front panel
2	Display	A 4 x 20 Liquid Crystal Display (LCD) shows the value of all current conditions including Temperature and Liquid Nitrogen Level. The display also shows any current alarm conditions that may exist.
3	Start Fill	Used to manually fill the freezer with liquid nitrogen. This button may also be pressed to clear “fog” from the storage area to increase visibility.
4	Stop Fill	Used to manually end the filling cycle. After pressing this button, the unit will not resume filling until the low-level setting is reached, or the “START FILL” button is pressed.
5	Setup	Used to adjust the value of all user adjustable functions including Level, Level Alarms, Temperature Alarms, Passwords, etc.
6	Alarm Mute	Used to silence the audible alarm. Also used to reset the latching alarm after the alarm condition is corrected.
7	Up Arrow	Used to increase number values during “Setup” modes. Press once to decrease incrementally. Hold the button to scroll through values quickly. May also be used to toggle “YES/NO” or “ON/OFF” values.
8	Down Arrow	Used to decrease number values during “Setup” modes. Press once to decrease incrementally. Hold the button to scroll through values quickly. May also be used to toggle “YES/NO” or “ON/OFF” values.
9	Escape	Used to exit any menu or setup function.
10	Enter	Used to Select any menu for editing or save any user adjustable setting.

TEC 3000 Quick Reference User Manual


Table 2: Bottom Panel Identification


NOTE: The recommended LN2 supply pressure is 22-35 psig (1.52 - 2.41 bar).

1	Temp A	Connection for Temperature A probe.
2	Temp B	Connection for Temperature B probe.
3	Power Connection	Connection for incoming power.
4	COM 1	Connection for RS-485 communications interface.
5	COM 2	Redundant Connection for RS-485 communications interface.
6	Manual Fill	Manual Fill button used to override the TEC 3000 software and manually fill the freezer. This button may also be pressed to clear "fog" from the storage area to increase visibility
7	Discreet Alarm Connections	These connections monitor the alarms for High Temp A, Low Level, Low Battery, and LN2 Supply. They can be wired to any alarm system that uses Normally Open (NO) contacts.
8	Global Alarm Connections	These connections monitor ALL alarms. They can be wired to any alarm system that uses Normally Open (NO) or Normally Closed (NC) contacts.
9	Main Harness Connection	Connection for the main wiring harness. All of the connections for the solenoid valves, lid switch, bypass sensor, etc. are located here.
10	Level Connection	Connection for level sensing system. Level sensing is achieved via differential pressure sensor, located inside the TEC 3000. The tube from the level sensing line on the freezer connects here. The connection is a 1/8" (3.2mm) OD hose barb.

Dewar Plumbing Connections

Connect a transfer line (included with freezer) from an LN₂ supply tank to the fill connection at the rear of the freezer. Optimum supply tank pressure is 22 to 35 psi (1.5 to 2.4 bar). Although the plumbing assembly has a 50 psi (3.45 bar) pressure relief device, it is recommended that the supply tank be pressurized below 35 psi (2.4 bar) to reduce the LN₂ “flash-off” rate during filling and to maximize the cryogenic valve life. The supply line can be insulated to minimize LN₂ transfer losses. After the transfer hose is securely coupled to the freezer and supply tank, ensure all connections are leak free by opening the valve of the LN₂ supply tank and apply a soap and water solution to each field joint. You should not see bubbles forming at any joint. Wipe away excess soap and water when finished. Before removing the transfer hose, ensure the LN₂ supply tank valve is closed. Slowly and carefully loosen the transfer hose connection to vent any remaining pressure in the line before disconnecting the hose.

NOTE: The recommended LN₂ supply pressure is 22-35 psig (1.52 - 2.41 bar).

TEC 3000 Quick Reference User Manual

NOTE: The recommended LN2 supply pressure is 22-35 psig (1.52 - 2.41 bar).

The following section describes how to adjust temperature alarm settings. At any time during the following procedure, the user may exit the level by pressing the “Escape” button repeatedly, until the display returns to the “monitor” display mode. After 30 seconds of inactivity, the controller will automatically return to the “monitor” display mode.

NOTE: Security Level 2 or higher is required to adjust temperature settings (see “Password and Security Setup” section for details).

To exit any menu screen and return to the previous menu, press “ESC” key.

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

High User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “ENTER”**

The display will read “Temperature Menus”.

Press Enter For
Temperature Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “Temp A Menu”.

NOTE: To adjust Temp B settings, press “Setup” instead of “Enter”.

Press Enter For
Temp A Menu
Or Press Setup For
Next Menu

4. **PRESS “ENTER”**

The display will read “Temp A Enabled”.

NOTE: Probe A may be disabled now by pressing the “▲ / ▼” keys. Press “Enter” to save settings.

Temp A
Enabled
Use ▲ / ▼ To Adjust
Press Enter To Save

5. **PRESS “SETUP”**

Until Display Reads “Temp A High Alarm”. Use the “▲ / ▼” keys to adjust the temperature. Hold either button to scroll quickly. Press “Enter” to save settings.

Temp A High Alarm
Use ▲ / ▼ To Adjust
Press Enter To Save
-110.0°C

6. **PRESS “SETUP”**

The display will read “Temp A Low Alarm”. Use the “▲ / ▼” keys to adjust the temperature. Hold either button to scroll quickly. Press “Enter” to save settings.

Temp A Low Alarm
Use ▲ / ▼ To Adjust
Press Enter To Save
-200.0°C

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

TEC 3000 Quick Reference User Manual

NOTE: The recommended LN2 supply pressure is 22-35 psig (1.52 - 2.41 bar).

The following section describes how to adjust liquid nitrogen level settings and the high/low level alarms.

NOTE: Security Level 2 or higher is required to adjust the Level and Level Alarm settings (See “Password and Security Setup” section for details).

NOTE: To exit any menu screen and return to the previous menu press “ESC” key.

HIGH LEVEL ALARM SETTING:

1. PRESS “SETUP”

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move the cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. PRESS “SETUP”

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

2. PRESS “ENTER”

The display will read “High Level Alarm”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save the settings. **NOTE:** No level settings can be adjusted to less than .5 inches (14 mm) within each other.

High Level Alarm
8.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

HIGH LEVEL SETTING:

1. PRESS “SETUP”

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. PRESS “SETUP”

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

3. PRESS “ENTER”

The display will read “High Level Alarm”.

High Level Alarm
8.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

4. PRESS “SETUP”

The display will read “High Level Setpoint”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

High Level Setpoint
7.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

TEC 3000 Quick Reference User Manual

Adjusting Level & Level Alarm Settings

NOTE: To exit any menu screen and return to the previous menu, press “ESC” key.

LOW LEVEL SETTING

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” arrow keys to scroll to the appropriate number. Press “Enter” to move the cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For Liquid
Level Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “High Level Alarm”.

High Level Alarm
8.0 In
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

Until the display reads “Low Level Setpoint”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

Low Level Setpoint
5.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

LOW LEVEL ALARM SETTING:

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Liquid Level Menus”.

Press Enter For
Liquid Level Menus
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “High Level Alarm”.

High Level Alarm
8.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

Until the display reads “Low Level Alarm”. Use the “▲ / ▼” keys to adjust the setting. Press “Enter” to save settings.

Low Level Alarm
4.0 in
Use ▲ / ▼ To Adjust
Press Enter To Save

TEC 3000 Quick Reference User Manual

Adjusting Display & Output Settings

The units of measurement displayed by the TEC 3000 may be adjusted to accommodate the needs of the user. Temperature measurement may be displayed in degrees Kelvin (°K), degrees Celsius (°C), or degrees Fahrenheit (°F). The amount of liquid nitrogen in the freezer may be displayed in inches (in), millimeters (mm). In addition, the amount of liquid nitrogen consumed by the freezer (liquid usage) may be shown on the display.

NOTE: Security Level 1 is required to adjust the display and output settings (See “Password and Security Setup”) section for details).

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.

Higher User Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “SETUP”**

Until the display reads “Display and Output”.

Press Enter For
Display And Output
Or Press Setup For
Next Menu

3. **PRESS “ENTER”**

The display will read “Temperature Units”. Use the “▲ / ▼” keys to toggle between °C (Celsius), °F (Fahrenheit), or °K (Kelvin). Press “Enter” to save settings.

Temperature Units
°C
Use ▲ / ▼ To Adjust
Press Enter To Save

4. **PRESS “SETUP”**

The display will read “Level Units”. Use the “▲ / ▼” keys to toggle between “in” (inches) or “mm” (millimeters). Press “Enter” to save settings.

Level Units
In
Use ▲ / ▼ To Adjust
Press Enter To Save

5. **PRESS “SETUP”**

The display will read “ Display Liquid Usage”. Use the “▲ / ▼” keys to select “ Yes” or “No”. Press “ Enter” to save settings.

After all adjustments are complete, wait 30 seconds and the TEC 3000 will return to the Main Display screen indicating new settings are in effect.

Calibration of Temperature Probes

Each of the temperature sensor probes (A and B) used with the TEC 3000 have been calibrated at the factory using the “Low Temperature Range” method. This calibration method provides a level of accuracy of +/-1.8°F (+/-1°C) when operated in the altitude range of 1000ft. to 1500ft. (305m to 457m). Further calibration should not be required unless desired by the end user. Refer to the TEC 3000 Technical Manual for information on calibration methods and procedures.

TEC 3000 Quick Reference User Manual

Password and Security Setup

The TEC 3000 can store up to 10 different passwords. Each password can be assigned its own security level ranging from Level 1 to Level 4. Table 3 below shows which settings can be changed with each security level. A security level of 4 is required to adjust any password. The default (or “Global”) password for the TEC 3000 is “3456”. All parameters may be adjusted by using this password. Record all passwords and security settings and store in a safe place. NOTE: MVE recommends changing the global password, as it is common to all units. If the global password has been forgotten, contact MVE Customer Service for details on how to reset passwords.

Table 3: Security Levels and Definitions

FEATURE	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4
Fill Start	X	X	X	X
Fill Stop	X	X	X	X
Alarm Mute	X	X	X	X
Change Display Units	X	X	X	X
Temp Settings		X	X	X
Level Settings		X	X	X
Time/Date		X	X	X
Calibration Probes		X	X	X
Change Languages		X	X	X
Hot Gas Bypass Settings		X	X	X
OFAF Setting			X	X
Communication Settings			X	X
Programming			X	X
Password Settings				X

1. **PRESS “SETUP”**

The controller will prompt for a password. Use the “▲ / ▼” keys to scroll to the appropriate number. If setting up passwords for the first time, use the default password of “3456”. Press “Enter” to move cursor to the next position in the password.

Higher Use Level
Required Use ▲ / ▼ To
Enter Password
0000

2. **PRESS “ENTER”**

The display will read “Temperature Menus”.

Press Enter For
Temperature Menus
Or Press Setup For
Next Menu

3. **PRESS “SETUP”**

Until the display reads “Password Menus”

Press Enter For
Password Menus
Or Press Setup For
Next Menu

4. **PRESS “ENTER”**

The display will read “Press Enter to change Global Password or press Setup for Next”.

Press Enter To
Change Global
Password Or Press
Setup For Next


NOTE: The default setting for the “Global” password is “3456”. The “Global” password may be used to change ALL settings for the TEC 3000.

Password and Security Setup

5. **PRESS “SETUP”**
To scroll to the desired password number.
6. **PRESS “ENTER”**
To select the desired password number.
7. **PRESS “▲” or “▼”**
To scroll to the appropriate number. Press “Enter” to move cursor to the next position in the password.
8. **PRESS “ENTER”**
After the desired password is selected. The display will now read “Password X Level”. Use the “▲ / ▼” keys to adjust the security level.
9. **PRESS “ENTER”**
Display will read “Confirm New Password?”. Use the “▲ / ▼” keys to select “YES” or “NO”. Selecting “YES” will save the new password and security setting. Selecting “NO” will abort all changes.
10. **PRESS “ENTER”**
To save settings.

Password 1
Use ▲ / ▼ To Adjust
Press Enter For Next
XXXX

Password X Level
Use ▲ / ▼ To Adjust
Press Enter For Next
Level X

Confirm New
Password?
NO
Use ▲ / ▼ To Adjust

At this time, the selected password and security level has been set and saved. Follow steps 7 - 11 to set additional passwords and security levels, or press “Escape” repeatedly to return to the display screen. After 5 minutes of inactivity, the controller will automatically return to the “monitor” mode. Use Table 4 to record passwords for future reference.

TEC 3000 Quick Reference User Manual

Alarms and Descriptions

NOTE: The recommended LN2 supply pressure is 22-35 psig (1.52 - 2.41 bar).

Table 4: Alarms and Descriptions

Alarm Display	Description
High Temp A	The temperature of Probe A is above the user defined High Temperature setting.
High Temp B	The temperature of Probe B is above the user defined High Temperature setting.
Low Temp A	The temperature of Probe A is below the user defined Low Temperature setting.
Low Temp B	The temperature of Probe B is below the user defined Low Temperature setting.
High Level	The depth of LN2 inside the freezer is above the user defined High level setting.
Low Level	The depth of LN2 inside the freezer is below the user defined Low level setting.
Usage Warning	The consumption of LN2 has doubled (See TEC 3000 Technical Reference Manual for details).
Usage Alarm	The consumption of LN2 has increased by a factor of 5 (See TEC 3000 Technical Reference Manual for details).
Fill Time	The amount of time required to complete a fill cycle exceeds the user defined Fill Time setting.
Bypass Time	The amount of time required to complete a bypass cycle exceeds the user defined Bypass Time setting.
Temp A Calibration	The temperature of Probe A is lower than absolute zero.
Temp B Calibration	The temperature of Probe B is lower than absolute zero.
Bypass Calibration	The temperature of the Bypass Probe is lower than absolute zero.
Low Battery	The voltage of the backup batteries has dropped below 21 volts.
Power Failure	The primary power has been disconnected for at least 60 minutes.
Lid Open	The lid on the freezer has been open longer than the user specified time.
Communication Loss	The controller has lost communications with the display.

If any alarms occur, contact your authorized MVE Distributor or customer / technical service.

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MVE Biological Solutions
 3055 Torrington Drive
 Ball Ground, GA 30107

MVE Inc. reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us as applicable. Copyright © 2019 MVE Inc., all rights reserved. MVE, MVE Fusion, QDRIVE, the MVE logo, and the MVE Security Through Systems logo are registered trademarks of MVE Inc. Other names may be trademarks of their respective owners.

MVE Biological Solutions
3055 Torrington Drive
Ball Ground, GA 30107