





MVE Vario™ TS Controller *(Touch Screen)*



Quick Reference Guide

MVE Vario TS Quick Reference Guide

Table of Contents

Product Identification

Display / Control Panel	2
Bottom Panel / Electrical / Physical Connections	3
Dewar Plumbing Connections.....	4
Adjusting Temperature Alarm Settings.....	5
Adjusting Inlet Temperature Settings	8
Adjusting Display and Output Settings.....	12
Password and Security Setup.....	14
Alarms and Definitions	17
Contact Information	17

Failure to follow Chart's best operating practices, as set forth in the operating manual, can result in loss of contents.

Active monitoring required.

Investigate all alarms or abnormal operating conditions immediately and address root causes.

Contact your Distributor for support.

Si no se siguen las prácticas operativas recomendadas de Chart, como se establece en el manual de operaciones, podría resultar en la pérdida de contenido.

Se requiere un control activo.

Investigue todas las alarmas o condiciones operativas anormales de inmediato y aborde las causas de origen.

Comuníquese con su distribuidor para obtener ayuda.

Il mancato rispetto delle migliori procedure operative di Chart, come descritto nel Manuale di funzionamento, può comportare la perdita di contenuto.

Richiesto monitoraggio attivo.

Esaminare immediatamente tutti gli allarmi o condizioni operative anomale e affrontare le cause radice.

Contattare il distributore per supporto.

Ne pas suivre les pratiques exemplaires d'exploitation de Chart, conformément au manuel d'utilisation, peut entraîner la perte de contenu.

Surveillance active nécessaire.

Enquêtez immédiatement sur toutes les alarmes ou les conditions d'exploitation anormales, et attaquez-vous aux causes profondes.

Contactez votre distributeur pour obtenir de l'assistance.

Die Nichtbeachtung der in dieser Betriebsanleitung aufgeführten bewährten Praktiken von Chart kann Inhaltsverluste nach sich ziehen.

Aktive Überwachung erforderlich.

Überprüfen Sie alle Alarme und ungewöhnlichen Betriebsbedingungen unverzüglich und beheben Sie die Ursache.

Wenden Sie sich bei Bedarf an Ihren Händler.

21323414 Rev A

MVE Vario TS Quick Reference Guide





Table 1: Front Panel Identification

Display	6" touchscreen, backlight
Freezer Status	Displays "IDLE", "BYPASSING", or "COOLING" based on the current freezer status
FOF CLEAR Key	Used to manually clear the fog from the storage area to increase visibility
END CLEAR Key	Used to manually terminate the fog clearing process
ALARM MUTE Key	Used to silence the audible alarm. Will reset the latching alarm once it has been corrected
SETUP Key	Used to access Setup Menus and parameters
Trend Graph	Visual X,Y graph parameter adjustable days and cooling temperature

MVE Vario TS Quick Reference Guide



1	Temp A Port	Connection for Temp A probe
2	Temp B Port	Connection for Temp B probe
3	Serial Number Barcode	Chart TS serial number written below barcode
4	30 VDC Power Input	Main power supply connection
5	Serial Port	RJ-45 connection for Serial/COM
6	Ethernet Port	Ethernet connection for networking
7	Global/Discrete Alarm Dry Contacts	Output connection for the remote monitoring of alarm conditions
8	Wire Harness Connection	12-pin wire harness connection to plumbing assembly, lid switch, and battery backup

MVE Vario TS Quick Reference Guide

Dewar Plumbing Connections

Connect a transfer line (included with freezer) from an LN2 supply tank to the fill connection at the rear of the freezer. Optimum supply tank pressure is 22 to 35 psi (1.5 to 2.4 bar). Although the plumbing assembly has a 50 psi (3.45 bar) pressure relief device, it is recommended that the supply tank be pressurized below 35 psi (2.4 bar) to reduce the LN2 “flash-off” rate during filling and to maximize the cryogenic valve life. The supply line can be insulated to minimize LN2 transfer losses. After the transfer hose is securely coupled to the freezer and supply tank, ensure all connections are leak free by opening the valve of the LN2 supply tank and apply a soap and water solution to each field joint. You should not see bubbles forming at any joint. Wipe away excess soap and water when finished. Before removing the transfer hose, ensure the LN2 supply tank valve is closed. Slowly and carefully loosen the transfer hose connection to vent any remaining pressure in the line before disconnecting the hose.



MVE Vario TS Quick Reference Guide

Adjusting Temperature A & B and Inlet Alarm Settings


The following section describes how to adjust temperature alarm settings. At any time during the following procedure, the user may exit the menu by pressing the “EXIT” button to return to the “monitor” display mode. After 60 seconds of inactivity, the controller will automatically return to the “monitor” display mode.

NOTE: Security Level 2 or higher is required to adjust temperature settings (see “Password and Security Setup” section for details).


To exit any menu screen and return to the previous menu, press “<” key.

1. Press “Setup”

Controller may prompt for a password. Type in the password using the number pad that appears and press “Enter”.



2. Press “Temperature Settings”



MVE Vario TS Quick Reference Guide

3. Press “Temperature A Settings”

NOTE: To access Temperature B Settings select “Temperature B Settings” instead.

CHART

TEMP A -85.9°C
TEMP B -81.1°C
CHAMBER -85.9°C
COOLING 17°C

IDLE

FOG CLEAR
END CLEAR
ALARM MUTE
SETUP

TEMPERATURE SETTINGS EXIT

-196.2°C

Temperature A Settings
Temperature B Settings
Inlet Temperature Settings
Cooling Temperature Settings

4. Press “ENABLED” or “DISABLED” next to “Temperature Probe A”

This will enable or disable the selected temperature probe. Pressing “ENABLED” will change the probe status to “DISABLED” and pressing “DISABLED” will change the probe status to “ENABLED”.

CHART

TEMP A -85.9°C
TEMP B -81.1°C
CHAMBER -85.9°C
COOLING 17°C

IDLE

FOG CLEAR
END CLEAR
ALARM MUTE
SETUP

TEMPERATURE A SETTINGS EXIT

Temperature Probe A ENABLED
High Alarm Setpoint -75.0°C
High Alarm ENABLED
Low Alarm Setpoint -95.0°C
Low Alarm ENABLED

Initiate High Temp. A Alarm Test Temperature A Calibration

MVE Vario TS Quick Reference Guide

5. Press the value displayed next to “High Alarm Setpoint”

The number pad will be displayed once the value to be adjusted is selected. Type in a new value for the High Alarm Setpoint using the number pad that appears on screen and then press “Enter” to save the new value. Be sure to include “-” when entering negative values.

CHART		TEMPERATURE A SETTINGS		EXIT	
TEMP A	-85.9°C	Temperature Probe A	ENABLED		
TEMP B	-81.1°C	High Alarm Setpoint	-75.0°C		
CHAMBER	-85.9°C	High Alarm	ENABLED		
COOLING	17°C	Low Alarm Setpoint	-95.0°C		
IDLE		Low Alarm	ENABLED		
FOG CLEAR		Initiate High Temp. A Alarm Test		Temperature A Calibration	
END CLEAR					
ALARM MUTE					
SETUP					

6. Press the value displayed next to “Low Alarm Setpoint”


The number pad will be displayed once the value to be adjusted is selected. Type in a new value for the Low Alarm Setpoint using the number pad that appears on screen and then press “Enter” to save the new value. Be sure to include “-” when entering negative values.

CHART		TEMPERATURE A SETTINGS		EXIT	
TEMP A	-85.9°C	Temperature Probe A	ENABLED		
TEMP B	-81.1°C	High Alarm Setpoint	-75.0°C		
CHAMBER	-85.9°C	High Alarm	ENABLED		
COOLING	17°C	Low Alarm Setpoint	-95.0°C		
IDLE		Low Alarm	ENABLED		
FOG CLEAR		Initiate High Temp. A Alarm Test		Temperature A Calibration	
END CLEAR					
ALARM MUTE					
SETUP					


MVE Vario TS Quick Reference Guide

Adjusting Inlet Temperature Settings

1. Press “Temperature Settings”



2. Press “Inlet Temperature Menus”



MVE Vario TS Quick Reference Guide

3. The current Hot Gas Bypass settings will be displayed

The current inlet temperature is displayed, along with the hot gas bypass setpoints.

CHART		INLET TEMPERATURE SETTINGS		EXIT
TEMP A	-85.9°C	Inlet Temperature	2.5°C	
TEMP B	-81.1°C	Hot Gas Bypass and Alarm	ENABLED	
CHAMBER	-85.9°C	Inlet Temperature Setpoint	-50.0°C	
COOLING	17°C	Supply Alarm Time Delay	5 minutes	
IDLE		Stuck Valve Alarm	DISABLED	
FOG CLEAR		Stuck Open Delay	30 minutes	Inlet Temperature Calibration
END CLEAR		Stuck Closed Delay	30 minutes	
ALARM MUTE				
SETUP				

4. Press “ENABLED” or “DISABLED” next to “Hot Gas Bypass and Alarm”

This will enable or disable the Hot Gas Bypass feature. Pressing “ENABLED” will change the hot gas bypass status to “DISABLED” and pressing “DISABLED” will change the hot gas bypass status to “ENABLED”.

CHART		INLET TEMPERATURE SETTINGS		EXIT
TEMP A	-85.9°C	Inlet Temperature	2.5°C	
TEMP B	-81.1°C	Hot Gas Bypass and Alarm	ENABLED	
CHAMBER	-85.9°C	Inlet Temperature Setpoint	-50.0°C	
COOLING	17°C	Supply Alarm Time Delay	5 minutes	
IDLE		Stuck Valve Alarm	DISABLED	
FOG CLEAR		Stuck Open Delay	30 minutes	Inlet Temperature Calibration
END CLEAR		Stuck Closed Delay	30 minutes	
ALARM MUTE				
SETUP				

MVE Vario TS Quick Reference Guide

5. Press the value displayed next to “Inlet Temperature Setpoint”

The number pad will be displayed once the value to be adjusted is selected. Type in a new value for the inlet temperature setpoint (if desired) using the number pad that appears on screen and then press “Enter” to save the new value.

INLET TEMPERATURE SETTINGS		EXIT
Inlet Temperature	2.5°C	
Hot Gas Bypass and Alarm	ENABLED	
Inlet Temperature Setpoint	-50.0°C	
Supply Alarm Time Delay	5 minutes	
Stuck Valve Alarm	DISABLED	
Stuck Open Delay	30 minutes	Inlet Temperature Calibration
Stuck Closed Delay	30 minutes	

6. Press the value displayed next to “Hot Gas Bypass Alarm Delay”

The number pad will be displayed once the value to be adjusted is selected. Type in a new value for the hot gas bypass alarm delay (if desired) using the number pad that appears on screen and then press “Enter” to save the new value.

INLET TEMPERATURE SETTINGS		EXIT
Inlet Temperature	2.5°C	
Hot Gas Bypass and Alarm	ENABLED	
Inlet Temperature Setpoint	-50.0°C	
Supply Alarm Time Delay	5 minutes	
Stuck Valve Alarm	DISABLED	
Stuck Open Delay	30 minutes	Inlet Temperature Calibration
Stuck Closed Delay	30 minutes	

MVE Vario TS Quick Reference Guide


Adjusting Display and Output Settings

The units of measurement displayed by the MVE TS may be adjusted to accommodate the needs of the user. Temperature measurement may be displayed in Kelvin (K), degrees Celsius (°C), or degrees Fahrenheit (°F).


NOTE: Security Level 1 is required to adjust the display and output settings (See “Password and Security Setup” section for details).

1. Press “Setup”

Controller may prompt for a password. Type in the password using the number pad that appears and press “Enter”.




2. Press “Display and Output Settings”



MVE Vario TS Quick Reference Guide

3. Press the units displayed next to “Temperature Units”

Press the “°C”, “°F”, or the “K” option for “Temperature Units”.



MVE Vario TS Quick Reference Guide

Password and Security Setup

The MVE TS can store up to 10 different passwords. Each password can be assigned its own security level ranging from Level 1 to Level 4. Table 3 below shows which settings can be changed with each security level. A security level of 4 is required to adjust any password. The default (or “Global”) password for the MVE TS is “3456”. All parameters may be adjusted by using this password. Record all passwords and security settings and store in a safe place. NOTE: MVE recommends changing the global password, as it is common to all units. If the global password has been forgotten, contact MVE Customer Service for details on how to reset passwords.

Table 3: Security Levels and Definitions

FEATURE	LEVEL 1	LEVEL 2	LEVEL 3	LEVEL 4
Fog Clear	X	X	X	X
End Clear	X	X	X	X
Alarm Mute	X	X	X	X
Change Display Units	X	X	X	X
Temp Settings		X	X	X
Time/Date		X	X	X
Calibration Probes		X	X	X
Change Languages		X	X	X
Hot Gas Bypass Settings		X	X	X
Communication Settings			X	X
Programming			X	X
Password Settings				X


MVE Vario TS Quick Reference Guide

This section details how to enable / disable password entry mode as well as how to change and setup multilevel security passwords.


NOTE: Security Level 4 is required to setup or change passwords.

1. Press “Setup”

Controller may prompt for a password. Type in the password using the number pad that appears and press “Enter”.



2. Press “Password Settings”



MVE Vario TS Quick Reference Guide

3. Press “ENABLED” or “DISABLED” next to “Password Entry Mode”

This will enable or disable the password entry mode. Pressing “ENABLED” will change the password entry mode to “DISABLED” and pressing “DISABLED” will change the password entry mode to “ENABLED”.

CHART		PASSWORD SETTINGS		EXIT
TEMP A	-85.9°C	Password Entry Mode		ENABLED
TEMP B	-81.1°C	Set Password 1	Set Password 6	
CHAMBER	-85.9°C	Set Password 2	Set Password 7	
COOLING	17°C	Set Password 3	Set Password 8	
IDLE		Set Password 4	Set Password 9	
FOG CLEAR		Set Password 5	Set Password 10	
END CLEAR		Set Global Password		
ALARM MUTE				
SETUP				

4. Press the desired password to setup or adjust.

CHART		PASSWORD SETTINGS		EXIT
TEMP A	-85.9°C	Password Entry Mode		ENABLED
TEMP B	-81.1°C	Set Password 1	Set Password 6	
CHAMBER	-85.9°C	Set Password 2	Set Password 7	
COOLING	17°C	Set Password 3	Set Password 8	
IDLE		Set Password 4	Set Password 9	
FOG CLEAR		Set Password 5	Set Password 10	
END CLEAR		Set Global Password		
ALARM MUTE				
SETUP				

5. Enter a new password and password level.

MVE Vario TS Quick Reference Guide

Alarms and Descriptions

Table 4: Alarms and Descriptions

Alarm Display	Description
High Temp A	The temperature of Probe A is above the user defined High Temperature setting.
High Temp B	The temperature of Probe B is above the user defined High Temperature setting.
Low Temp A	The temperature of Probe A is below the user defined Low Temperature setting.
Low Temp B	The temperature of Probe B is below the user defined Low Temperature setting.
Cooling Time	The amount of time required to complete a cool cycle exceeds the user defined Cool Time setting.
Bypass Time	The amount of time required to complete a bypass cycle exceeds the user defined Bypass Time setting.
Temp A Calibration	The temperature of Probe A is lower than absolute zero.
Temp B Calibration	The temperature of Probe B is lower than absolute zero.
Bypass Calibration	The temperature of the Bypass Probe is lower than absolute zero.
Low Battery	The voltage of the back up batteries has dropped below 21 volts.
Power Failure	The primary power has been disconnected for at least 60 minutes.
Lid Open	The lid on the freezer has been open longer than the user specified time.
Communication Loss	The controller has lost communications with the display.



**If any alarms occur, contact your authorized
MVE Distributor or customer / technical service.**

Customer/Technical Service:

USA:

Phone: 1-844-683-2796

Fax: 1-470-552-2200

Europe:

Phone: +44 (0) 1344 403 100

Fax: +44 (0) 1344 429 224

Asia:

Phone: +61 297 494333

Fax: +61 297 494666



MVE Biological Solutions
3055 Torrington Drive
Ball Ground, GA 30107
+1-844-683-2796
USA www.mvebio.com

**Intended Use & Indication for Use
for Cryogenic Storage and/or Transport**

STORAGE ONLY

MVE Freezers are intended for the maintenance of cryogenic temperatures during storage for the indication of preserving human or animal biological products, samples, or specimens (e.g., blood, blood products, cells, tissues, etc.) during storage.

STORAGE AND TRANSPORT

MVE Dewars and Vapor Shippers are intended for the maintenance of cryogenic temperatures during storage or transportation for the indication of preserving human or animal biological products, samples, or specimen: (e.g., blood, blood products, cells, tissues, etc.) during storage and or transportation.



Ref 21124996 Rev D