



OPERATION MANUAL

Condition Monitoring Systems CryoBeacon & SmartTag





© 2025 MVE Biological Solutions reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Register your logger:
<https://register.mvebio.com>

Customer Service:
+1-844-683-2796
monitor@mvebio.com

**Visit www.mvebio.com
for more information.**

Table of Contents

○ Table of Contents	3
GENERAL	4
○ Purpose	4
○ Distributor Identification	4
Products Covered	4
Regulatory Information	4
○ Warranty	4
CRYOBEACON	5
○ Basic Safety Instructions	5
○ Functional Description	5
Design and Components	5
○ Operation	6
○ Tips for Saving Energy	6
○ Cleaning & Care	6
○ Service Instructions	6
Connecting External Probe	6
Probe Adjustment or Replacement	7
○ Technical Data	11
○ Housing Dimensions & Mounting Template	12
○ Storage Intervals	12
SMARTTAG	13
○ Basic Safety Instructions	13
○ Functional Description	13
Design & Components	13
Display Description	14
○ Operation	14
Startup	14
Start Recording	15
Stop Recording	15
Export PDF Report	15
Checking Battery Status & Charging	15
Tracking Functions	16
Data Recording & Transmission Interval	16
○ Cleaning & Care	16
○ Service Instructions	17
Probe Adjustment or Replacement	17
Charging the Logger	21
○ Technical Data	23
○ Orientation	24
○ Shock Values (G Values Related to Drop in Height)	24
○ Shock Sensors	25
○ Battery Life with Transmission Intervals	25
TROUBLESHOOTING	26
CERTIFICATE OF CONFORMANCE	27

GENERAL

Purpose

This operating manual contains important information for safe, correct and economical operation of the MVE Condition Monitoring Systems (CMS), CryoBeacon and SmartTag. Observing the instructions will help to avoid dangers, reduce repair and increase reliability and prolong the service life of the loggers.

Distributor Identification

MVE Biological Solutions US, LLC
3055 Torrington Drive
Ball Ground, Georgia 30107
USA

Telephone: +844-683-2796
Email: monitor@mvebio.com
Website: www.mvebio.com

Products Covered

This manual applies to the following products:

- CryoBeacon – Passive condition monitoring for shipment and storage tracking.
- SmartTag – Near real-time condition monitoring for shipment tracking (temperature, light, shock, orientation, acceleration, and geo-position).

Regulatory Information

These products comply with applicable international standards and directives.

CryoBeacon

RED Directive 2014/53/EU
ROHS Directive 2011/65/EU and WEEE Directive 2012/19/EU
REACH Directive 2009/1907/EU
CE Mark Low Voltage 2014/35/EU
FCC CFR Title 47 Part 15 Subpart B
Contains FCC ID: X8WBT840
ICES-003 Contains IC: 4100A-BT840

SmartTag

FCC CFR Title 47 Part 15 Subpart B
CE Marking: These devices conform to the low voltage directive (2014/35/EU)
Compliant with RTCA/DO-160G
UN38.3
ROHS Directive 2011/65/EU and WEEE Directive 2012/19/EU



Warranty

MVE will warrant condition monitoring systems (CMS) to be free from defects in material or workmanship for one (1) year from initial shipment to the original purchaser from the MVE warehouse.

MVE's liability for goods that do not conform to the product warranty shall be strictly limited, at MVE's sole option, to one of the following:

- (i) Repair of the goods;
- (ii) Replacement of the goods with new goods that conform to the product warranty; or
- (iii) Refund of the price paid by the original buyer for the goods.

Repair, replacement, or refund, at MVE's sole option, shall be the original buyer's sole and exclusive remedy for any actual failure of the goods to conform to the product warranty. MVE shall have no liability to the original buyer for any loss, damage, cost, or expense of any nature whatsoever, except as expressly provided herein.

Any CMS replaced under warranty shall be warranted for the greater of: (a) the unexpired portion of the original warranty period applicable to the device, or (b) ninety (90) days from the date of replacement.

For warranty claims, contact Customer/Technical Support at +844-683-2796 or monitor@mvebio.com.

CRYOBEACON

Basic Safety Instructions

Never open the housing of the logger. Defective parts or components must only be replaced by parts and components approved by MVE Biological Solutions.

If the logger shows visible damage, it must not be put into operation.

The RJ45 connector only works with proprietary probes provided or approved by MVE Biological Solutions and is not supported by any standard.

Functional Description

The CryoBeacon logs ambient temperature and humidity. With the external analog probe, cryogenic temperature conditions are monitored during transit or in storage. The data can be accessed via the MVECloud after scanning with the Tec4App.

Design and Components



1	Logger Serial Number
2	QR Code for Data Readout
3	Temperature & Humidity Sensor
4	Mounting Holes
5	RJ-45 Jack for External Probe



Operation

Before operating, ensure the unit is dry and the battery charge is >10% (approximately 3 weeks of battery charge remaining). The battery charge for a specific logger can be viewed on the app or website after scanning the logger.

Do not place the logger near open flames or other heat sources (strong sunlight, gas oven, heater, etc.).

Ensure the temperature/humidity sensor and membrane are not covered.

Do not immerse the logger in liquid.

Avoid throwing or skidding the logger or applying excessive external forces.

If no external probe is connected, cover the RJ45 port to keep it dry and protected from liquids.

Operate, store and distribute the logger in the following ambient conditions only: min. -29°C/-20°F; max. +60°C/+140°F.

If CMS is being installed on a charged shipper or cold dewar, NER should not be measured until 5 hours after installation to allow the cork/cover to equilibrate with the dewar which was properly cooled per MVE's Vapor shipper and Dewar instructions.

Tips for Saving Energy

For optimal battery life, initiate Bluetooth data readouts 1-2 times per day.

The operation of the logger in low temperatures may decrease the battery runtime.

Cleaning & Care

Before cleaning, disconnect the external probe from the logger, if applicable. Ensure the connector port is kept dry and protected from all liquids. Clean only the exterior surface of the logger using a soft cloth that is lightly dampened with water. Do not immerse the logger in water, expose it to running water, or allow moisture to enter the device. Do not use harsh chemicals, solvents, abrasive cleaning agents, or sharp objects, as these may damage the logger and void the warranty.

Service Instructions

Connecting External Probe

Remove the RJ45 protection cover, if applicable.

Plug the external probe into the RJ port.

Pull gently on the probe to ensure proper connection.

Note: The locking tab of the connector must face upwards.

After connecting the external Pt-100 temperature probe, the CryoBeacon records both internal temperature and temperature of the probe. If no probe is connected, the CryoBeacon only records the internal temperature.

Probe Adjustment or Replacement

Tools & Materials Required

#1 Phillips Screwdriver

GE Silicone Sealant: Advanced Silicone, Clear, 3 oz, Tube, 301% to 500% Elongation Range

1. Loosen the (4) quarter-turn fasteners using a #1 Phillips screwdriver by turning the fastener 90 degrees counterclockwise to remove the cover of the CMS.




2. Remove the probe by pressing the tab on the cable connection and pulling outward, away from the datalogger.



3. Using a #1 or #2 Phillips screwdriver, remove the two SS #6-32 x 3/16" LG 18-8 screws securing the CryoBeacon to its bracket. Always remove the CryoBeacon before repositioning the probe to avoid probe damage and failed temperature readings.




4. Remove existing silicone from the probe/cork mating surfaces by pulling it off gently.



5. Remove the wrapped probe cable from the wire retention tab & remove the cable tie if needed.



6. To remove the probe, push the exposed probe tip into the cork and then pull the cable from the top side of the assembly. See the direction in the image below




- 7.** To adjust the height of the probe, see instructions below:

Disclaimer: The standard probe protrusion recommended by MVE Biological Solutions from the bottom of the cork foam to the tip of the probe:


- SC 4/2V: 1 5/8" (41.3mm)
- SC 4/3V: 1 1/2" (38.1mm)

It is the responsibility of the end user to independently verify and validate any dimensions to ensure suitability for their intended application.

- a.** Pull downward on the probe from the cork if the desired height is to be lower in the dewar or vapor shipper.



- b.** Push upward on the probe into the cork, while also pulling on the cable in the cover if the desired height is to be higher in the dewar or vapor shipper.



- 8.** With the probe in the desired location, place the wrapped cable back underneath the wire retention clip.



- 9.** Replace the CryoBeacon onto the bracket as shown below. Replace the screws (2x, circled in red). Tighten by hand. If using an alternate method of tightening, the screws should not exceed a maximum torque of 8 in-lbs (0.9 Nm). For alignment of the CryoBeacon, the edge of the CryoBeacon case and the edge of the bracket should be flush (indicated by the green arrows). This will create the proper angle to match the window on the cover of the housing.



- 10.** Plug the wire connector back into the CryoBeacon as shown below.



- 11.** Replace the cover of the housing back by aligning the locating tab on the cover with the corresponding slot on the lower housing. See below in the red circle.



- 12.** Tighten the quarter turns with a screwdriver by turning each 90 degrees clockwise to secure the cover of CMS. If using an alternate method of tightening, the screws should not exceed a maximum torque of 8 in-lbs (0.9 Nm).



- 13.** Apply silicone sealant around the probe (or wire) where it meets the cork to seal the gap between the probe or wire and the hole in the cork. Apply the silicone 360 degrees around the probe, and it is recommended not to have any strands of silicone protruding from the silicone bead. After applying, allow 24 hours for the silicone to fully cure before putting into service.




Technical Data


	Description	Details
Record & transmission configuration	Record interval	• 10 minutes
	Transmission interval	• On-demand when scanning the CryoBeacon
Operating conditions	Operating temperature	• -29°C to +60°C
	Storage temperature	• +5°C to +40°C (recommended)
	Humidity range	• 0% rH to 100% rH, non-condensing
Liquid Nitrogen (LN₂)	Down to -200°C	<ul style="list-style-type: none"> • Analog Pt-100 4-wire class A probe • Diameter: Cable 2.1mm; Sheath 3mm • Length: Cable 1000mm; Sheath 100mm • RJ45 connector • Temp. range: -200°C ... +100°C
Ambient Sensors calibrated by manufacturer	Temperature	<ul style="list-style-type: none"> • -40°C to +85°C with accuracy of ±0.2°C • 0.1°C resolution
	Humidity	<ul style="list-style-type: none"> • 0% rH ... 100% rH non-condensing • ±2% rH accuracy (at 20% rH ... 80% rH) • 0.1% rH resolution
Memory & Logging	Non-volatile memory event and interval triggered with thresholds	• Capacity: 5,000 data points in circular buffers
Data transfer & analysis	Wireless	• Data transfer via 2.4 GHz Bluetooth connection
2.4 GHz for data transfer	2.4 GHz to read out data	• Max. 400m range (line of sight)
Battery	CR2477 3V Lithium coin cell	<ul style="list-style-type: none"> • Battery life depends on the data rate and ambient temperature; approximately 12 months • For transport, no marking obligation of the lithium metal button cell (ADR special provision 188 and IATA DGR)
Housing & Mounting	ABS Housing Mounting: e.g. Screw M3 ISO 7380 FL or industrial adhesive tape	<ul style="list-style-type: none"> • Dimensions: 88mm x 47mm x 19mm • Distance of mounting holes: 80mm • Maximum tightening torque: 8 in-lbs (0.9 Nm) • Weight: approx. 47g
Water Ingress	IP 50 protection	<ul style="list-style-type: none"> • Sensor protection through membrane • Penetration of fluids is to be prevented (Corrosion damage/short circuit)
Approvals & Standards	EC Declaration of Conformity (CE marking)	<ul style="list-style-type: none"> • RED Directive 2014/53/EU • ROHS Directive 2011/65/EU and WEEE • REACH Directive 2009/1907/EU • CE Mark Low Voltage 2014/35/EU
	FCC	<ul style="list-style-type: none"> • FCC CFR Title 47 Part 15 Subpart B • Contains FCC ID: X8WBT840
	IC (Canada)	• ICES-003 Contains IC: 4100A-BT840
Export information	Product tariff code	• 9031 8080
	Country of origin	• DE (EU)
	Data regarding legal control	• ECCN = N, AI = no

Housing Dimensions & Mounting Template

Housing Cross Section
(dimensions in mm)



Housing Measurements (dimensions in mm)
Note: Print this page as a mounting template



Storage Intervals

Record Interval (min)	Full memory reached in (days)
1	~4 days
2	~8 days
3	~12 days
4	~15 days
5	~19 days
6	~23 days
10	~38 days
12	~46 days
15	~57 days
20	~76 days
30	~114 days
60	~229 days

SMARTTAG

Basic Safety Instructions

Never open the housing of the logger. Defective parts or components must only be replaced with parts and components approved by MVE Biological Solutions.

If the logger shows visible damage, it must not be put into operation.

Disconnect the logger from the power supply before any cleaning and maintenance and after each use.

Only charge using the MVE Biological Solutions USB charging cable.

Only connect the logger with the charging cable to a 5V USB power supply.


Functional Description

The SmartTag is a multifunctional and flexible real-time tracking solution. Integrated sensors measure temperature, light, shock, orientation, acceleration, and geo-position and send them to the cloud software via the mobile network.


The SmartTag can be used as a stand-alone near real-time monitoring system working as a live logger.

The SmartTag is a multi-use logger and can be recharged and used multiple times. It operates as the link to the MVECloud, offering full IoT functionality via mobile connectivity worldwide.

Design & Components



Display Description



Operation

Before use, make sure the charging cable and probe connector are dry. Only charge the logger at ambient temperature of 0°C to +45°C.

If the logger shows visible damage, do not put into operation..

Do not operate the logger near open flames or other sources of heat.

Protect the logger and cables from heat and moisture.

Do not immerse the logger in liquid.

Avoid throwing or skidding the logger or applying excessive external forces.

Operate, store and distribute the logger in the following ambient conditions only: min. -29°C/-20°F; max. +60°C/+140°F.


If CMS is being installed on a charged shipper or cold dewar, NER should not be measured until 5 hours after installation to allow the cork/cover to equilibrate with the dewar.


Startup

The SmartTag is deactivated prior to shipping. A short press of either the STOP or START button will illuminate the LCD display. The SmartTag will display the temperature on the LCD screen, but it is not recording data at this point.



Start Recording

Press and hold the START button for at least 5 seconds to initiate the record mode. The record icon  will illuminate on the left side of the LCD display. The recorded data can be viewed in the MVECloud or Tec4app.




Stop Recording

The STOP button is deactivated on the device to prevent accidental stoppage during shipment.

Recording can be set up to stop automatically with checkpoints/destinations or manually in the cloud.

After recording is stopped, MAX and MIN temperature can be viewed on the device by pressing the START button to scroll through the data.



NOTE: USB connection to a computer or charging the logger will also stop data recording.

Export PDF Report

The PDF reports can be obtained any time via the cloud which can be accessed via any web browser. Optionally, the SmartTag can be connected to a computer via the Micro USB port at the bottom of the logger to retrieve the PDF date report.

NOTE: USB connection to a computer or charging the logger will also stop data recording.

Checking Battery Status & Charging


A short press will illuminate the LCD screen to view the battery icon in the upper right corner of the screen.

If the battery power is low before the first use, fully charge the logger using the provided charging cable.

While charging, the flash icon will be displayed next to the battery icon.

See “Charging” under “Service Instructions” in this manual.

NOTE: USB connection to a computer or charging the logger will also stop data recording.



Approved Airlines

Email Customer Service for list of approved airlines.

Data Recording & Transmission Interval

The SmartTag measures temperature, shock, light and geo-position.

The SmartTag records the measured parameters in an adjustable interval. The default interval is 10 minutes.

The SmartTag transmits the parameters to the cloud in an adjustable interval. The default interval is 60 minutes.

The intervals can be changed by the MVE on request (shorter intervals = higher power consumption, longer intervals = lower power consumption).

Cleaning & Care

Clean only the exterior surface of the logger using a soft cloth that is lightly dampened with water. Do not immerse the logger in water, expose it to running water, or allow moisture to enter the device. Do not use harsh chemicals, solvents, abrasive cleaning agents, or sharp objects, as these may damage the logger and void the warranty.

Service Instructions

Probe Adjustment or Replacement

Tools & Materials Required


#1 Phillips Screwdriver

GE Silicone Sealant: Advanced Silicone, Clear, 3 oz, Tube, 301% to 500% Elongation Range

1. Loosen the (4) quarter-turn fasteners using a #1 Phillips screwdriver by turning the fastener 90 degrees counterclockwise to remove the cover of CMS.
2. Remove probe by turning the lock nut counter-clockwise on the cable connection to loosen, then pull the two ends of the connection to disconnect.




3. Lift the SmartTag out from its bracket by pulling either end of the datalogger. Always remove the SmartTag before repositioning the probe to avoid probe damage and failed temperature readings.
4. Remove existing silicone from the probe/cork mating surfaces by pulling it off gently.



- 5.** Remove the screws (2 ea SS #6-32 X 3/16" LG 18-8) that mount the bracket to the housing and remove the bracket. This will allow the probe wire to move freely.



- 6.** To remove the probe, push the exposed probe tip into the cork and then pull the cable from the top side of the assembly. See the direction in the image below





- 7.** To adjust the height of the probe, see instructions below:

Disclaimer: The standard probe protrusion recommended by MVE Biological Solutions from the bottom of the cork foam to the tip of the probe:

- SC 4/2V: 1 5/8" (41.3mm)
- SC 4/3V: 1 1/2" (38.1mm)

It is the responsibility of the end user to independently verify and validate any dimensions to ensure suitability for their intended application.

- a.** Pull downward on the probe from the cork if the desired height is to be lower in the dewar or vapor shipper.
- b.** Push upward on the probe into to cork, while also pulling on the cable in the cabinet if the desired height is to be higher in the dewar or vapor shipper.



- 8.** With the probe in the desired location, fold the wire, replace the cable tie, and place the wrapped portion of the cable where the wire retention clip will hold the cable in place. Route the connection side of the wire following the red line shown below.



- 9.** Replace the bracket into the housing, by aligning the location holes in the bracket with the locating pins on the housing. The screw holes in the bracket should line up with the threaded inserts in the housing. Once aligned, place the screws (2x) through the bracket into the threaded inserts. Tighten by hand. If using an alternate method of tightening, the screws should not exceed a maximum torque of 8 in-lbs (0.9 Nm).



- 10.** Place the SmartTag back into the bracket, with the display & buttons facing upward. Route the cable as shown below in red.



- 11.** Reconnect the probe by aligning the arrows on the connector(s). This connection is keyed so no excessive force should be applied when making the connection, as it may damage the pins in the connector. Tighten the nut to secure the connection.



- 12.** Replace the cover of the housing back by aligning the locating tab on the cover with the corresponding slot on the lower housing. See below in the red circle.



- 13.** Tighten the quarter turns with a screwdriver by turning each 90 degrees clockwise to secure the cover of CMS. If using an alternate method of tightening, the screws should not exceed a maximum torque of 8 in-lbs (0.9 Nm).




- 14.** After the CMS is fully assembled, the SmartTag must be restarted to be used in shipment and data recording. See SmartTag “Operation” for instructions to “Start Recording” to restart data recording in this manual.


Charging the Logger

1. Loosen the (4) quarter-turn fasteners using a #1 Phillips screwdriver by turning the fastener 90 degrees counterclockwise to remove the cover of CMS.
2. Pull the white tab on the bottom of the SmartTag to remove the dust cover of the charging port.



3. Plug the charging cable into the SmartTag. The screen will flash "PDF" indicating that the logger is searching for a computer to upload a PDF document via USB. When using a standard power supply and not a computer, this can be ignored and the screen will switch back to normal screen where it will display a  symbol, indicating that the unit is charging. A red LED will also appear, showing that the unit is charging.



4. When fully charged, the screen will display a  symbol indicating that the unit is fully charged. The red LED will turn off, and a blue LED will turn on, also indicating that the unit is fully charged. When fully charged, remove the cable and replace the dust cover on the charging port.



5. Replace the cover of the housing back by aligning the locating tab on the cover with the corresponding slot on the lower housing. See below in the red circle.




6. Tighten the quarter turns with a screwdriver by turning each 90 degrees clockwise to secure the cover of CMS. If using an alternate method of tightening, the screws should not exceed a maximum torque of 8 in-lbs (0.9 Nm).



Technical Data

	Description	Details
Record & Transmission	Record interval	• 1 / 2 / 5 / 10 / 30 / 60 min
	Transmission interval	• 10 / 20 / 30 / 60 / 120 min
Accelerometer	3 axes: x, y and z 1-8 g	<ul style="list-style-type: none"> • 1-8 g factory-calibrated • Resolution: 1 g
Operating conditions	Operating temperature	• -29°C to +60°C
	Storage temperature	• +5°C to +40°C (recommended)
	Humidity range	• 0% rH to 100% rH, non-condensing
Temperature Sensor	Calibrated by manufacturer	<ul style="list-style-type: none"> • -30°C to 70°C • ±0,5°C (-20°C to +40°C) • 0.1 °C resolution
Light Sensor	Detect whether box is open	• 0 to 52.000 Lux
Geo-Positioning	Detect geolocation of device	• Geo-positioning via cellular and WiFi triangulation
Battery	Li-Ion	• Battery runtimes on following pages
Casing	ABS casing	• Dimensions: 103 x 61.3 x 33 mm
Weight	SmartTag	• 300 g
GSM Module Config	Global roaming included	<ul style="list-style-type: none"> • Global LTE 4G/5G • Fallback 2G network
Approvals & Standards	Compliance	<ul style="list-style-type: none"> • FCC CFR Title 47 Part 15 Subpart B • CE Marking: These devices conform to the low voltage directive (2014/35/EU) • RTCA/DO-160G compliant • UN38.3 • ROHS Directive 2011/65/EU and WEEE Directive 2012/19/EU
Protection	IP 64	• Sensor protection through membrane
Versions	Pt100 sensor 4-wire connection, class A	• 100 mm husk & 1000 mm cable

Orientation




Shock Values (G Values Related to Drop in Height)

Sensor Threshold	10cm	30cm	50cm	70cm	90cm	110cm	130cm	150cm
Lev1	✓	✓	✓	✓	✓	✓	✓	✓
Lev2	✗	✓	✓	✓	✓	✓	✓	✓
Lev3	✗	✗	✓	✓	✓	✓	✓	✓
Lev4	✗	✗	✗	✓	✓	✓	✓	✓
Lev5	✗	✗	✗	✗	✓	✓	✓	✓
Lev6	✗	✗	✗	✗	✗	✓	✓	✓
Lev7	✗	✗	✗	✗	✗	✗	✓	✓
Lev8	✗	✗	✗	✗	✗	✗	✗	✓

Shock Sensors

The shock sensor indicates the impact and the time when the impact occurred.



Battery Life with Transmission Intervals

Battery runtime in days is equal to the transmission interval in minutes.

Transmission Interval (min)	Battery Runtime (days)
120	120
60	60
30	30
20	20
10	10

TROUBLESHOOTING

Logger	Fault	Resolution
CryoBeacon	App not recognizing the logger	<ul style="list-style-type: none"> Verify smart device has Bluetooth enabled.
	Abnormal temperature readings	<ul style="list-style-type: none"> Visually inspect the external probe and/or wire for damage.
	Temperature not as expected inside the dewar	<ul style="list-style-type: none"> Verify the probe length to ensure probe is in the desired location.
	Temperature reading inside the dewar is reading ambient	<ul style="list-style-type: none"> Ensure the connection of the external probe to the CryoBeacon is securely installed.
	CMS cover is loose	<ul style="list-style-type: none"> Confirm all quarter-turn fasteners are secure.
SmartTag	Data is not recording to the cloud	<ul style="list-style-type: none"> Verify the R appears on the LCD screen. See "Start Recoding" instructions in this manual.
	Abnormal temperature readings	<ul style="list-style-type: none"> Visually inspect the external probe and/or wire for damage.
	Temperature not as expected inside the dewar	<ul style="list-style-type: none"> Verify the probe length to ensure probe is in the desired location.
	CMS cover is loose	<ul style="list-style-type: none"> Confirm all quarter-turn fasteners are secure.
	Light sensor is not sensing light	<ul style="list-style-type: none"> Remove the cover and verify the logger bracket is secure in the housing.

Tutorials are available at www.mvebio.com

Contact Customer Service at 844-683-2796 or email monitor@mvebio.com for further assistance.

CERTIFICATE OF CONFORMANCE

 <div style="display: inline-block; vertical-align: middle; margin-left: 10px;"> <h3 style="margin: 0;">Certificate of Conformance</h3> </div>			
<p>This document certifies that the accompanied Conditional Monitoring System or component meets MVE Biological Solutions' internal quality standards, is calibrated, and operates as designed.</p> <p>MVE Biological Solutions has been audited and certified to EN ISO 13485: 2016</p>			
Sign:		Date:	Sep 9th, 2025
Name:	Max Mortensen, Ph.D.		
Title:	Vice President Quality Assurance & Regulatory Affairs		
<p>This certificate is not to be interpreted as a warranty, neither expressed nor implied.</p>			

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference in radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. The Federal Communications Commission (FCC) warns the users that changes or modifications to the unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Americas

Customer Service:
+1-844-683-2796
monitor@mvebio.com

EMEA

Customer Service:
+44 7718 488236
monitor@mvebio.com

China (Domestic)

Customer Service:
+86 28 6572 9660
Jessy Yang
Direct: +86 28 6572 9668
monitor@mvebio.com


Asia Pac Rim

Customer Service:
+01 404-383-1900
monitor@mvebio.com

**Visit www.mvebio.com
for more information.**

Global Network, Local Impact: With MVE's global footprint, we have the reach to support customers worldwide. Yet, our approach remains personal, with local expertise and support teams dedicated to your success.

Our Locations



MVE Biological Solutions
3055 Torrington Drive
Ball Ground, GA 30107
Ph +1-470-552-2500
Toll Free +1-844-683-2796
Fax +1-404-383-1901
monitor@mvebio.com